


Special fluorescent lamps


Medical/therapy



Low-pressure mercury discharge lamps producing (U.V.) radiation for medical purposes. The lamp is mainly used to treat hyperbilirubinaemia (/52), and is only available for medical and industrial use.

Applications


- Treatment of neonatal jaundice (hyperbilirubinaemia)
- Climatologic testing equipment



Dimensions in mm

Type	Cap/base	Lumen output lm	Ordering number
/52			
'TL' 20W /52	G13	360	8222 206 80700





Spectrum of a blue lamp TL/52: emission peak at wavelength of 450 nm.